

Anti-Ihh Rabbit antibody



Catalog Number: 167661

Product name

Anti-Ihh Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Ihh Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Ihh

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/2000

IHC: 1/20

Validations



Western blot detection of Ihh in
HepG2,U2OS,HL-60,Lncap using Ihh Rabbit
mAb(1:1000 diluted)