

15 Lipoxygenase 2 antibody

Catalog Number: 107661

Product name

15 Lipoxygenase 2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

15 Lipoxygenase 2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in MCF7 cells, HEK-293 cells. Observed molecular weight by Western-blot: 76kd

Preparation

This antibody was obtained by immunization of 15 Lipoxygenase 2 recombinant protein (Accession Number: NM_001141). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107661 (ALOX15B antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:107661 (ALOX15B antibody) at dilution of 1:100 (under 10x lens)