

Anti-ARF3 / ADP-ribosylation factor 3 antibody



Catalog Number: 102361

Product name

Rabbit IgG

Anti-ARF3 / ADP-ribosylation factor 3 antibody

Applications

Immunogen

ELISA, WB, IF, ICC/IF

[Human ARF3 / ADP-ribosylation factor 3 \(His Tag\) recombinant protein](#)

Dilutions

Specificity

WB: 5-10 µg/ml

Human ARF3 / ADP-ribosylation factor 3

ELISA: 0.1-0.2 µg/ml

Antibody description

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human ARF3 / ADP-ribosylation factor 3.


Rabbit polyclonal to ARF3 / ADP-ribosylation factor 3

ICC/IF: 0.5-1.5 µg/mL

Preparation

Validations

Produced in rabbits immunized with purified, recombinant Human ARF3 / ADP-ribosylation factor 3 (rh ARF3 / ADP-ribosylation factor 3; P61204; Met1-Lys181). ARF3 / ADP-ribosylation factor 3 specific IgG was purified by Human ARF3 / ADP-ribosylation factor 3 affinity chromatography.



Lanes	A
Items	
Sample (whole cell lysate)	MCF7
Sample Volume (µg/lane)	30
Gel	13%SDS-PAGE reducing gel
Recommended Concentration	5-10 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L) at 1:5000 dilution. <small>Developed using Odyssey imaging system.</small>
Explanation	Predicted band size : 20 kDa Observed band size : 17 kDa

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.


Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

ARF3 / ADP-ribosylation factor 3 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



ARF3 / ADP-ribosylation factor 3 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence

Anti-ARF3 / ADP-ribosylation factor 3 antibody



Catalog Number: 102361

Immunofluorescence staining of ARF3 in HeLa cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-human ARF3 polyclonal antibody (1 µg/ml) at 4°C

overnight. Then cells were stained with the Alexa Fluor®594-conjugated Goat Anti-rabbit IgG secondary antibody (red) and counterstained with DAPI (blue). Positive staining was localized to cytoplasm.