

NUDT3 antibody



Catalog Number: 113361

Product name

NUDT3 antibody

Specificity

human, mouse; other species not tested.

Antibody description

NUDT3 Rabbit Polyclonal antibody. Positive WB detected in mouse cerebellum tissue, mouse liver tissue. Observed molecular weight by Western-blot: 20-30kd

Preparation

This antibody was obtained by immunization of NUDT3 recombinant protein (Accession Number: NM_006703). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:113361(NUDT3 antibody) at dilution of 1:800