

# Anti-PRSS3 antibody



Catalog Number: 176361

## Product name

Polyclonal

Anti-PRSS3 antibody

## Ig Type

## Specificity

Rabbit IgG

Human

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to PRSS3

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human trypsin-3. aa range:117-166

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality



Fig1:: Western Blot analysis of various cells using Trypsin-3 Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000