

Barkor-Specific antibody

Catalog Number: 108361

Product name

Barkor-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Barkor-Specific Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in mouse brain tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_014924). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IF: 1:10-1:100


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108361(Barkor-Specific antibody) at dilution of



IP Result of anti-Barkor-Specific (IP:Catalog No:108361, 4ug; Detection:Catalog No:108361 1:500) with HeLa cells lysate 2800ug.



Immunofluorescent analysis of MCF-7 cells using Catalog No:108361(Barkor-Specific Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated

Barkor-Specific antibody



Catalog Number: 108361

AffiniPure Goat Anti-Rabbit IgG(H+L)