

E.coli rppH/nudH/ygdP/b2830/JW2798 recombinant protein



Catalog Number: 518361

General Information

Gene Name Synonym

RNA pyrophosphohydrolase, (Di)nucleoside polyphosphate hydrolase, Ap5A pyrophosphatase, rppH

Protein Construction

A DNA sequence encoding the E.coli rppH/nudH/ygdP/b2830/JW2798 Met1-Gly176 is expressed.

Organism

E.coli

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 50mM Tris, 500mM NaCl, 10% glycerol, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



E.coli rppH/nudH/ygdP/b2830/JW2798 recombinant protein