

Human PAPSS1/ATPSK1/PAPSS (His tag) recombinant protein



Catalog Number: 519361

General Information

Gene Name Synonym

Bifunctional 3'-phosphoadenosine 5'-phosphosulfate synthase 1, PAPS synthase 1, PAPSS 1, Sulfurylase kinase 1, SK 1, ATPSK1, PAPSS

Protein Construction

A DNA sequence encoding the Human PAPSS1/ATPSK1/PAPSS Met1-Ala624 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM Tris-HCl buffer (pH 8.0), 100mM NaCl and 20% glycerol.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PAPSS1/ATPSK1/PAPSS (His tag)
recombinant protein