

Anti-ZNF473 antibody



Catalog Number: 175761

Product name

Anti-ZNF473 antibody

Specificity

Human, Mouse, Rat, Dog, Cow, Horse, Rabbit

Antibody description

Rabbit polyclonal antibody to ZNF473

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human zfp100/znf473 701-830/871

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF

Dilutions

WB:1:500-2000

IHC-P:1:400-800

IHC-F:1:400-800

IF:1:100-500

ICC/IF:1:100-500

Validations




Fig1: Sample: Bone (mouse) Lysate at 40 ug; Primary: Anti- NF473 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 100 kD; Observed band size: 100 kD