

# Anti-hnRNP U (1B11) Mouse antibody



Catalog Number: 168071

## Product name

Anti-hnRNP U (1B11) Mouse antibody

## Specificity

Human

## Antibody description

hnRNP U (1B11) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human hnRNP U protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of hnRNP U in K562, Jurkat and HeLa cell lysates and using hnRNP U mouse mAb(1:1000 diluted). Predicted band size: 110KDa. Observed band size: 110KDa.