

# FAM129B antibody



Catalog Number: 110471

## Product name

FAM129B antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FAM129B Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of FAM129B recombinant protein (Accession Number: BC067366). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110471(FAM129B Antibody) at dilution of 1:1000

Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:110471(FAM129B Antibody) at dilution of 1:50 (under 40x lens)