

Anti-CD33 antibody



Catalog Number: 142471

Product name

Anti-CD33 antibody

Immunogen

[Human CD33 / Siglec-3 \(His Tag\) recombinant protein](#)

Specificity

Human CD33 / Siglec-3

Antibody description

Rabbit Monoclonal to Human CD33

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CD33 / Siglec-3 and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



CD33 / Siglec-3 Antibody (APC), Rabbit MAb, Flow Cytometry