

# Anti-MYC antibody



Catalog Number: 104471

## Product name

Anti-MYC antibody

## Specificity

Recognize N-terminal in fusion proteins.

## Antibody description

Mouse monoclonal to MYC

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, a synthetic Myc tag peptide conjugated to KLH. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1



## Applications

WB

## Dilutions

WB: 1/5000-1/10000

## Validations



Anti-Myc Tag Antibody (N-terminal Myc), Mouse MAbs, Western blot