

SNX15 antibody



Catalog Number: 115471

Product name

SNX15 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SNX15 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, COLO 320 cells, human testis tissue. Observed molecular weight by Western-blot: 50kd and 38kd

Preparation

This antibody was obtained by immunization of SNX15 recombinant protein (Accession Number: NM_013306). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:115471(SNX15 antibody) at dilution of 1:400