

PHOX2B antibody



Catalog Number: 107471

Product name

PHOX2B antibody

Immunogen

[Human PHOX2B Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

PHOX2B Mouse Monoclonal antibody. Positive WB detected in Neuro-2a cells. Positive IHC detected in neuroblastoma tissue. Observed molecular weight by Western-blot: 32-36 kDa

Preparation

This antibody was obtained by immunization of PHOX2B recombinant protein (Accession Number: NM_003924). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:107471(PHOX2B Antibody) at dilution of



1:2500

Immunohistochemistry of paraffin-embedded neuroblastoma tissue slide using (PHOX2B Antibody) at dilution of 1:500 (under 40x lens).