

# Anti-ICAM1 (4F5) Mouse antibody



Catalog Number: 168471

## Product name

Anti-ICAM1 (4F5) Mouse antibody

## Specificity

Human

## Antibody description

ICAM1 (4F5) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CD54(ICAM-1) protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


WB, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/200

## Validations



Immunocytochemistry staining of SW480 cells fixed by anhydrous methanol for 2 h at -20°C and using anti-CD54(ICAM-1) mouse mAb



(dilution 1:200).

Western blot detection of CD54(ICAM-1) in SW480 and Raji cell lysates using CD54(ICAM-1) mouse mAb (1:1000 diluted). Predicted band size:96KDa.Observed band size:96KDa.