

COQ6 antibody

Catalog Number: 109471

Product name

COQ6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

COQ6 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 51kd

Preparation

This antibody was obtained by immunization of COQ6 recombinant protein (Accession Number: NM_182476). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109471(COQ6 antibody) at dilution of 1:200