

# CYP1B1 Polyclonal Antibody



Catalog Number: 161871

---

## Product name

CYP1B1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to CYP1B1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 244-543 of human CYP1B1 (NP\_000095.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - CYP1B1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CYP1B1 antibody at dilution of 1:100 (40x lens).  
Blue: DAPI for nuclear staining.