

# Anti-SLC25A12 Rabbit antibody



Catalog Number: 168871

## Product name

Anti-SLC25A12 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, IP

## Antibody description

SLC25A12 Rabbit monoclonal antibody

## Dilutions

1:1000

1:20

## Preparation

Antigen: A synthetic peptide of human SLC25A12

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of SLC25A12 in Rat Brain, C6, 3T3 cell lysates using SLC25A12 Rabbit mAb (1:1000 diluted). Predicted band size: 75 kDa. Observed band size: 75 kDa.