

Anti-SLC25A12 Rabbit antibody



Catalog Number: 168871

Product name

Anti-SLC25A12 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IP

Antibody description

SLC25A12 Rabbit monoclonal antibody

Dilutions

1:1000

1:20


Preparation

Antigen: A synthetic peptide of human SLC25A12

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of SLC25A12 in Rat Brain,C6,3T3 cell lysates using SLC25A12 Rabbit mAb(1:1000 diluted).Predicted band size:75kDa.Observed band size:75kDa.

Clonality

Monoclonal