

TDG Polyclonal Antibody



Catalog Number: 161171

Product name

TDG Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to TDG

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TDG (NP_003202.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - TDG Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using TDG antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.