

Mouse CRISP-1 / CRISP1 (His Tag) recombinant protein



Catalog Number: 502171

General Information

Gene Name Synonym

Acidic epididymal glycoprotein 1; Sperm-coating glycoprotein 1

Protein Construction

A DNA sequence encoding the mouse Crisp1 (NP_033768.3) (Met1-His244) was expressed with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 20

Molecular Mass

The recombinant mouse Crisp1 consists 236 amino acids and predicts a molecular mass of 27.1 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CRISP-1 / CRISP1 Protein (His Tag) SDS-PAGE