

# LCK antibody



Catalog Number: 112171

---

## Product name

Polyclonal

LCK antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

LCK Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, Jurkat cells. Observed molecular weight by Western-blot: 50-60 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of LCK recombinant protein (Accession Number: BC013200). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:112171(LCK Antibody) at dilution of 1:1000

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality