

PRCP antibody

Catalog Number: 114171

Product name

PRCP antibody

IP: 1:200-1:1000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PRCP Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human liver tissue. Positive IP detected in L02 cells. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 56 kDa

Validations




Preparation

This antibody was obtained by immunization of PRCP recombinant protein (Accession Number: NM_005040). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114171 (PRCP antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114171 (PRCP antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000




Immunohistochemical of paraffin-embedded

PRCP antibody

Catalog Number: 114171

human lung cancer using Catalog
No:114171(PRCP antibody) at dilution of 1:50

IP Result of anti-PRCP (IP:Catalog No:114171,
4ug; Detection:Catalog No:114171 1:300) with
L02 cells lysate 1800ug.



(under 40x lens)