

CDON antibody



Catalog Number: 109171

Product name

CDON antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CDON Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells. Observed molecular weight by Western-blot: 140kd

Preparation

This antibody was obtained by immunization of CDON recombinant protein (Accession Number: XM_011542863). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:109171(CDON antibody) at dilution of 1:1000