

RCAN2 antibody



Catalog Number: 114571

Product name

RCAN2 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RCAN2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HT-1080 cells, mouse brain tissue. Positive IHC detected in human gliomas tissue, human skin tissue. Observed molecular weight by Western-blot: 24kd, 35kd

Preparation

This antibody was obtained by immunization of RCAN2 recombinant protein (Accession Number: BC038509). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114571(RCAN2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114571(RCAN2 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114571(RCAN2 antibody) at dilution of 1:100 (under 40x lens)