

UNC119 antibody



Catalog Number: 116571

Product name

UNC119 antibody

Specificity

Human, Mouse, Rat, Dog; other species not tested.

Antibody description

UNC119 Rabbit Polyclonal antibody. Positive IF detected in MDCK cells.

Preparation

This antibody was obtained by immunization of UNC119 recombinant protein (Accession Number: NM_005148). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of MDCK cells using Catalog No:116571(UNC119 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)