

Anti-FUBP1/FBP Rabbit antibody



Catalog Number: 168571

Product name

Anti-FUBP1/FBP Rabbit antibody

Applications

WB, IHC-P, ICC/IF, FC, IP

Specificity

Human, Mouse, Rat

Dilutions

1:1000

Antibody description

FUBP1/FBP Rabbit monoclonal antibody

1:500

Preparation

Antigen: A synthetic peptide of human FUBP1

1:50

1:50

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

1:20

Validations

Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG



Western blot detection of FUBP1/FBP in K562, Rat Brain, HeLa cell lysates using FUBP1/FBP Rabbit mAb(1:1000 diluted). Predicted band size: 68kDa. Observed band size: 74kDa.