

Human CD153/CD30L/TNFSF8 recombinant protein



Catalog Number: 500971

General Information

Gene Name Synonym

CD30 ligand

Protein Construction

A DNA sequence encoding the human CD30L (NP_001235.1) (Gln63-Asp234) was expressed with two additional amino acids (Gly & Pro) at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 63

Molecular Mass

The recombinant human CD30L consists of 172 amino acids and predicts a molecular mass of 19.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CD153 / CD30L / TNFSF8 Protein