

Human BCL2L11 / Bim (His Tag) recombinant protein



Catalog Number: 501971

General Information

Gene Name Synonym

Apoptosis regulator Bcl-X

Protein Construction

A DNA sequence encoding the human BCL2L11 (NP_612815.1) (Ala2-Arg120) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human BCL2L11 consists of 130 amino acids and predicts a molecular mass of 15

KDa. It migrates as an approximately 16-20 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50mM Tris, 10% Glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human BCL2L11 / Bim Protein (His Tag) SDS-PAGE