

LAG3 Polyclonal Antibody



Catalog Number: 163971

Product name

LAG3 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to LAG3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 29-170 of human LAG3 (NP_002277.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:1000 - 1:4000

Validations



Western blot - LAG3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using LAG3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.