

Human IL1RL2/IL1RRP2 (Fc Tag) recombinant protein



Catalog Number: 503971

General Information

Gene Name Synonym

IL-36 receptor; Interleukin-1 receptor-related protein 2

Protein Construction

A DNA sequence encoding the human IL1RL2 (NP_003845.2) (Met1-Tyr337) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 20

Molecular Mass

The recombinant human IL1RL2 consists 559 amino acids and predicts a molecular mass of 63.2 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL1RL2 / IL-1Rrp2 Protein (Fc Tag)