

Revdofilimab reference antibody



Catalog Number: 176971

Product name

Revdofilimab reference antibody

Specificity

Human TNFRSF4 / OX40 / CD134

Antibody description

Anti-TNFRSF4 / OX40 / CD134 Antibody (Revdofilimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-TNFRSF4 / OX40 / CD134 Antibody (revdofilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFRSF4 / OX40 / CD134 Antibody (revdofilimab) is more than 95%, determined by SEC-HPLC.