

ZWINT antibody



Catalog Number: 117271

Product name

ZWINT antibody

Immunogen

[Human ZWINT Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZWINT Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in human skin tissue. Observed molecular weight by Western-blot: 31kd

Preparation

This antibody was obtained by immunization of ZWINT recombinant protein (Accession Number: NM_007057). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



human skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:117271(ZWINT antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using ZWINT antibody Catalog No:117271 at 1:100 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent DNA dye).