

Mouse CD24 (Fc Tag) recombinant protein



Catalog Number: 501671

General Information

Gene Name Synonym

Cell surface retention sequence-binding protein 1;
Extracellular link domain-containing protein 1

Protein Construction

A DNA sequence encoding the mouse CD24 (Q8BHC0) (Met1-Arg52) was expressed, fused with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asn 27

Molecular Mass

The recombinant mouse CD24 /Fc is a disulfide-

linked homodimer. The reduced monomer comprises 267 amino acids and has a predicted molecular mass of 29.7 kDa. The apparent molecular mass of the protein is approximately 48-52 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD24 / Ly-52 / CD24A Protein (Fc Tag)
SDS-PAGE