

# Human IL20RA/IL-20RA (Fc Tag) recombinant protein



Catalog Number: 503371

## General Information

### Gene Name Synonym

Cytokine receptor class-II member 8; Cytokine receptor family 2 member 8; IL-20R1; ZcytoR7

### Protein Construction

A DNA sequence encoding the human IL20RA (NP\_055247.3) (Met1-Lys250) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Val 30

### Molecular Mass

The recombinant human IL20RA consists 459 amino acids and predicts a molecular mass of 52 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human IL20RA / IL-20RA Protein (Fc Tag) SDS-PAGE