

Anti-VN1R2 antibody



Catalog Number: 176371

Product name

Anti-VN1R2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to VN1R2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human vn1r2. aa range:88-137

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1::; Western Blot analysis of various cells using V1RL2 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000