

# Anti-PDK2 Mouse antibody



Catalog Number: 168371

## Product name

Anti-PDK2 Mouse antibody

WB: 1/500 - 1/2000

FC: 1/200 - 1/400

## Specificity

Human

ELISA: 1/10000

## Antibody description

PDK2 Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human PDK2 (AA: 178-404) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type


Mouse IgG2b

## Applications


WB, FC

## Dilutions

## Validations



Flow cytometric analysis of A431 cells using PDK2 mouse mAb (green) and negative control (red).



Western blot detection of PDK2 in A549, MCF7, C6 and HeLa cell lysates using PDK2 mouse mAb (1:2000 diluted). Predicted band size: 46 kDa. Observed band size: 46 kDa.