

GRIA3 Polyclonal Antibody



Catalog Number: 160771

Product name

GRIA3 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GRIA3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 200-490 of human GRIA3 (NP_000819.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC, IF


Dilutions

Validations



Western blot - GRIA3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GRIA3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - GRIA3 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using GRIA3 antibody . Blue: DAPI for nuclear staining.