

COX8A Polyclonal Antibody



Catalog Number: 164081

Product name

COX8A Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to COX8A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-69 of human COX8A (NP_004065.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:1000 - 1:2000

Validations



Western blot - COX8A Polyclonal Antibody

Western blot analysis of extracts of mouse skeletal muscle, using COX8A antibody .