

Anti-HtrA2/Omi (8G11) Mouse antibody



Catalog Number: 168081

Product name

Anti-HtrA2/Omi (8G11) Mouse antibody

Specificity

Human

Antibody description

HtrA2/Omi (8G11) Mouse monoclonal antibody

Preparation

Antigen: Recombinant human HtrA2/Omi protein.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB, IP

Dilutions


WB: 1/500

Validations



Western blot detection of HtrA2/Omi in 293T cell lysates using HtrA2/Omi mouse mAb(dilution 1:500). Predicted band size:49kDa. Observed band

size:36kDa.



Immunoprecipitation analysis of HeLa cell lysates using HtrA2/Omi mouse mAb.