

Anti-WWP1 Rabbit antibody



Catalog Number: 169081

Product name

WB: 1/500-1/2000

Anti-WWP1 Rabbit antibody

IHC: 1/25-1/100

Specificity

Human

Antibody description

WWP1 Rabbit polyclonal antibody

Preparation

Antigen: Fusion protein of human WWP1

Formulation

-20°C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC-P

Dilutions

Validations



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane: PC3 cell,
Primary antibody: 169081(WWP1 Antibody) at
dilution 1/750, Secondary antibody: Goat anti
rabbit IgG at 1/8000 dilution, Exposure time: 10



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 169081(WWP1 Antibody) at dilution 1/20. (Original magnification: ×200)