

# Anti-WWP1 Rabbit antibody



Catalog Number: 169081

## Product name

Anti-WWP1 Rabbit antibody

WB: 1/500-1/2000

IHC: 1/25-1/100

## Specificity

Human

## Antibody description

WWP1 Rabbit polyclonal antibody

## Preparation

Antigen: Fusion protein of human WWP1

## Formulation

-20°C, pH7.4 PBS, 0.05% NaN<sub>3</sub>, 40% Glycerol

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


ELISA, WB, IHC-P

## Dilutions

## Validations



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane: PC3 cell,  
Primary antibody: 169081(WWP1 Antibody) at  
dilution 1/750, Secondary antibody: Goat anti  
rabbit IgG at 1/8000 dilution, Exposure time: 10



The image is immunohistochemistry of paraffin-  
embedded Human esophagus cancer tissue using  
169081(WWP1 Antibody) at dilution 1/20.  
(Original magnification: ×200)