

CCT2 antibody

Catalog Number: 109081

Product name

CCT2 antibody

Immunogen

[Human CCT2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

CCT2 Rabbit Polyclonal antibody. Positive IF detected in NIH/3T3 cells. Positive WB detected in mouse testis tissue, Jurkat cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 50-58 kDa

Preparation

This antibody was obtained by immunization of CCT2 recombinant protein (Accession Number: NM_006431). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109081(CCT2 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed NIH/3T3 cells using Catalog No:109081(CCT2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)