

SNX32 antibody



Catalog Number: 115481

Product name

SNX32 antibody

Immunogen

[Human SNX32 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

SNX32 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of SNX32 recombinant protein (Accession Number: NM_152760). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115481(SNX32 Antibody) at dilution of 1:600