

Phospho-PRKACA-T197 pAb



Catalog Number: 166481

Product name

Phospho-PRKACA-T197 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-PRKACA-T197

Preparation

Antigen: A synthetic phosphorylated peptide around T197 of human PRKACA (NP_002721.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-PRKACA-T197 pAb

Western blot analysis of extracts of 293T cells, using Phospho-PRKACA-T197 antibody at 1:1000 dilution. 293T cells were treated by Insulin (100nM) for 10 minutes or treated by PMA/TPA (200nM) for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.