

Anti-BNIP3 Rabbit antibody



Catalog Number: 167481

Product name

Anti-BNIP3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

BNIP3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human BNIP3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

ICC/IF: 1/50

Validations



Immunocytochemistry of BNIP3(green) in HeLa cells using BNIP3 Rabbit pAb at dilution 1/50, and



Western blot detection of BNIP3 in Rat Brain lysates using BNIP3 Rabbit pAb(1:500 diluted). Predicted band size:28kDa. Observed band size:28kDa.