

Anti-SCGN / Secretagogin antibody



Catalog Number: 100881

Product name

Anti-SCGN / Secretagogin antibody

Immunogen

[Human SCGN / Secretagogin \(His Tag\) recombinant protein](#)

Specificity

Human SCGN / Secretagogin

Antibody description

Mouse monoclonal to SCGN / Secretagogin

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human SCGN / Secretagogin (rh SCGN / Secretagogin; NP_008929.2; Asp2-Pro276). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

IHC-P

Dilutions

IHC-P: 10-30 µg/mL

Validations



SCGN / Secretagogin Antibody, Mouse MAb, Immunohistochemistry

Immunochemical staining of human SCGN in human pancreas with mouse monoclonal antibody (15 µg/mL, formalin-fixed paraffin embedded sections).