

GATA3 antibody



Catalog Number: 110881

Product name

GATA3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GATA3 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells. Positive IHC detected in human breast cancer tissue, human colon cancer tissue. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of GATA3 recombinant protein (Accession Number: NM_002051). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110881(GATA3 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:110881(GATA3 Antibody) at dilution of 1:50 (under 40x lens)