

JAKMIP2 antibody



Catalog Number: 111881

Product name

Polyclonal

JAKMIP2 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

JAKMIP2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 95kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of JAKMIP2 recombinant protein (Accession Number: NM_014790). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111881(JAKMIP2 antibody) at dilution of 1:1500

Clonality